

=====  
Scope: The scope of this document is to highlight changes made for TRAC2ES Mobile version 1.0.7.1.

This release contains enhancements & fixes of the following change requests (CRs):

Known Issues:

Installing TRAC2ES MOBILE

In some cases, the install progress bar may go away with the appearance that nothing is happening, or the install has failed. Please, be patient, the install is still running in the background.

Installing TRAC2ES MOBILE

When using the DII COE Installer, the setup TRAC2ES MOBILE.exe is not properly called.

To complete a DII COE install, the setup TRAC2ES MOBILE.exe has to be invoked manually.

Enhancements (since v1.0.6.0):

CR689 – Removed the installation of Netscape 4.7.x from the TRAC2ES MOBILE installer and updated any text referencing Netscape 4.7.x.

CR533 – ICD9 codes have been updated from the database with this release. The ability to update ICD9 through the reference table module has also been added.

CR674 – Commercial ticket number and cost have been added to the movement tab of the PMR as well as being added to the attendant tab.

CR676 – Allergy information has been moved from the Vitals tab to the movement tab.

CR677 – The rate field on the Oxygen tab has been expanded to accept up to 30 characters.

CR678 – Added a new List box to the IV tab called the IV Location Type and having the values Heparin Lock, Central, Peripheral-Central Line.

CR679 – Updated the rules for the Litter/Ambulatory movement types.

CR680 – Added SP02 Rate to the Vital tab and added SP02 rate comments to the oxygen tab.

CR681 – Added the following fields to the new Arterial Blood Gases section of the Vitals tab (pH, P02, PC02, HC03, BE).

CR682 – Added the following fields to the new Ventilator Settings section of the Oxygen tab (Mode, FI02, TV, Rate, PEEP, PS).

Bug Fixes (since v1.0.6.0):

CR539 – Updated Reference Table to eliminate updating each reference table twice as well as removing the duplicate notifications regarding the update.

CR627 – Fixed the problem where TRAC2ES MOBILE was loading the PMR's seq\_pmr\_id into the Person ID field.  
CR665 – The crew type has been fixed so the text properly displays correctly. CR666 – Only allowing users to edit PMRs in the submitted or re-submitted states.

#### Enhancements:

CR50 - Allow the drop lists to behave as they do in Trac2es Web. Base the list selection

on the keystroke of the user. CR51 - Set ITV locations similar to how it's done in TWeb. CR437 - Give the user the ability to change their user account information without

having to shut down TRAC2ES MOBILE.

#### Bug Fixes:

CR151 - No error received for having more than 50 characters in Unit fields.

CR152 - Missing visual indicators of format requirements. Visual indicators on

not displayed for several fields on the PMR Clinical tab (Body Temperature, Pulse, Blood Pressure, Weight, Glucose, Pulse Oximeter Setting, and Oxygen Delivery Rate). CR153 - Rate field on IV tab is not validating field length constraint. CR154 - No check on fields requiring a limit to the number of characters entered. CR160 - Missing columns in Mission Reply tables. Patient table is missing personnel status. Attendant table is missing ID Type and Personnel Status. CR187 - Sorting medical specialties by Bed Report dates fail on MTF screen. CR195 - Nationality saved as "N/A" is returned as default (United States) on edit. CR444 - Casualty Event search with date range is not inclusive. CR535 - Space Type is always set to ambulatory when converting a mass PMR to a full PMR. CR537 - When a mass PMR is retrieved for editing theater settings are not returned. CR538 - Missing cite number is PMR header information. CR540 - Not closing key synchronization message

boxes on initial display causes issue.

Enhancements (since v1.0.4.3):

- CR59 - Provide the capability for a TRAC2ES MOBILE user to set/reset their password. The "Change Password" tab has been added to TRAC2ES MOBILE Administration Screen, which allows a user to enter a new password.
- CR305 - Install TRAC2ES MOBILE using InstallShield. TRAC2ES MOBILE is now installed by InstallShield. Significant changes have been made to the Readme.txt as a result.
- CR438 - Automate the key synchronization process. TRAC2ES MOBILE will automatically synchronize keys when needed. A user no longer needs to perform this function manually. The Security Administration Tab has been removed from the TRAC2ES MOBILE Administration Screen. A user can still invoke the manual process by clicking on the "Execute Key Request" button.
- CR549 - Upgrade browser packaged with TRAC2ES MOBILE to Netscape 4.7.8. Netscape 4.7.8 is the browser now packaged with TRAC2ES MOBILE. Other Browsers will work also now but didn't when Trac2es 1<sup>st</sup> came on line.
- CR552 - Provide the capability import/export PMRs using TRAC2ES MOBILE. The "Save Draft" icon is now used for exporting PMRs, A user can now export a PMR by clicking on the export icon located on toolbar of PMR screens. A user can now import a PMR through the "File" menu, "Import Message Request" menu-item.

Bug Fixes (since v1.0.4.3):

- CR165 - Issues viewing reports when TRAC2ES MOBILE is installed on D:\ A user can select the browser they wish to use for reports via the TRAC2ES MOBILE Administration Screen, "Preferences" tab.
- CR550 - Certain video cards cause systems to lock up on TRAC2ES MOBILE shut-down. If this problem persists, refer to CR550 Alternative Instructions.
- CR551 - Some special characters cause reports to not be returned.
- CR155 - Report Date/Time is not returned on TRAC2ES MOBILE reports.

CR550 Alternative Instructions:

Due to interoperability issues that exist between the Java Platform and certain video cards, some Mobile clients have experienced system lock ups and general system degradation upon shut down of TRAC2ES MOBILE.

A list of problematic video cards is below. While this list is non-inclusive of all problematic cards, it does include those that have been identified.

- ATI Radeon 7500
- ATI Radeon 7200 VE
- ATI Rage 128 GLAGP
- ATI Radeon 7000
- ATI Mobile M6
- ATI Mobility FireGL 7800
- ATI Radeon Mobility 7500
- ATI RAGE MOBILITY AGP
- ATI RAGE MOBILITY-M1 AGP2X

RAGE MOBILITY M7  
Nvidia GeForce 256 AGP  
ATI Radeon 9700

Your video card can be identified by conducting the following steps:

For Windows 2000 users:

- 1 From the start menu, navigate to "Settings" and select "Control Panel".
- 2 Double-click "System".
3. Click the "Device Manager" button. Page 3  
RELEASE\_NOTES.txt
- 3 From the Device Manager Console, expand any lists necessary until a "Display Adapters" list is present.
- 4 Expand the "Display Adapters" list to identify your video card.

For Windows NT users:

- 1 From the start menu, navigate to "Settings" and select "Control Panel".
- 2 Double-click "Display".
- 3 Select the "Settings" tab.
- 4 Click the "Display Type" button, your video card information will be in the "Adapter Type" section.

If your computer contains a listed video card or you experience any of the above symptoms, you will need to launch the TRAC2ES MOBILE client from the alternative startup file, "run TIRUIS Alternate.bat", located in h\tracmc\bin

For additional information regarding TRAC2ES MOBILE, see README.txt

